


CONSERVE PPE

Setup

View the PPE videos at med.emory.edu/PPE



After gown and glove removal...



1 **Sanitize** hands.

2 **Put on** fresh gloves.

3 **Place** wipe on table.

5 **Remove** eyewear and **place** on wipe.

6 **Sanitize** gloves.

7 **Remove** procedure mask, lower strap and then upper. **Store** mask.

8 **Sanitize** gloves.

Face Shield

9 **Wipe** front and back of shield.

10 **Wipe** elastic band.

11 **Wipe** foam band.

12 **Wipe** table.

13 **Place** shield upside down to dry.

14 **Sanitize** gloves.

15 **Remove** gloves.

16 **Wash hands** with soap and water.

Goggles

9 **Wipe** front and back of lens.

10 **Wipe** both ear pieces.

11 **Wipe** table.

12 **Place** goggles on clean table to dry.

13 **Sanitize** gloves.

14 **Remove** gloves.

15 **Wash hands** with soap and water.